




**FEATHER SAFETY RAZOR CO., LTD. (Japan)**

| PRODUCT CODE | PRODUCT DESCRIPTION                                                                                                                                                                                     |                                                                                      |
|--------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------|
| M10/8060002  | <b>Feather Microtome Blade No. R-35</b> (50 pieces/box)<br>Material: Stainless Steel<br>Dimension (T)x(W)x(L) : 0.254 x 8 x 80 Edge Angle 35"<br>Edge Finish: Resin Coating Platinum Spattering         |   |
| M10/8010002  | <b>Feather Microtome Blade No. S-35</b> (50 pieces/box)<br>Material: Stainless Steel<br>Dimension (T)x(W)x(L) : 0.254 x 8 x 80 Edge Angle 35"<br>Edge Finish: Resin Coating Platinum Spattering         |   |
| M10/8007002  | <b>Feather Microtome Blade "High Profile"</b> (50 pieces/box)<br>Material: Stainless Steel<br>Dimension (T)x(W)x(L) : 0.31 x 14 x 75.7 Edge Angle 35"<br>Edge Finish: Resin Coating Platinum Spattering |  |
| M09/8200002  | <b>Feather Microtome Holder No. 130A with Blade (No. S-35)</b><br>(per set)                                                                                                                             |                                                                                      |