






FEATHER SAFETY RAZOR COMPANY (Japan)

PRODUCT CODE	PRODUCT DESCRIPTION	
0FE-M10/8060002	Feather Microtome Blade No. R-35 (50 pieces/box) Material: Stainless Steel Dimension (T)x(W)x(L) : 0.254 x 8 x 80 Edge Angle 35" Edge Finish: Resin Coating Platinum Spattering	
0FE-M10/8010002	Feather Microtome Blade No. S-35 (50 pieces/box) Material: Stainless Steel Dimension (T)x(W)x(L) : 0.254 x 8 x 80 Edge Angle 35" Edge Finish: Resin Coating Platinum Spattering	
0FE-M10/8007002	Feather Microtome Blade "High Profile" (50 pieces/box) Material: Stainless Steel Dimension (T)x(W)x(L) : 0.31 x 14 x 75.7 Edge Angle 35" Edge Finish: Resin Coating Platinum Spattering	
0FE-M09/8200002	Feather Microtome Holder No. 130A with Blade (No. S-35) (per set)	