



TRILUX
NEUES LICHT.

TRILUX GmbH CO. & KG (Germany)

PRODUCT DESCRIPTION

Trilux LED OR-Light - Aurinio L 120

For ceiling installation with 1 luminary body 120.000 lux.
Consisting of: luminary body, gimbal, spring arm, cantilever boom arm, flange tube and ceiling flange for rooms with height 3.000mm until max. 3.500mm.



Trilux LED OR-Light - Aurinio L 160

For ceiling installation with 1 luminary body 160.000 lux.
Consisting of: luminary body, gimbal, spring arm, cantilever boom arm, flange tube and ceiling flange for rooms with height 3.000mm until max. 3.500mm.

